



# DNA methylation: Illumina 27K Array


Ali Balubaid, Azari Bantan, Turki Sobahy



# Workflow Overview



# BeadChip Technology



# Bisulfite Conversion

CpG island

----- ATTCTAG**CG**


----- TAAGAT**CG**

If methylated

----- ATTCTAG**CG**

----- TAAGAT**CG**

Bisulfite cytosine deamination




If **NOT** methylated

----- ATTCTAG**U**


----- TAAGAT**CG**

# Fragment Hybridization



If methylated




If **NOT** methylated



# Fragment Extension



# Allele Specificity



# Illumina's Infinium HumanMethylation27 BeadChip

27,578 CpG dinucleotides

14,495 genes

1  $\mu$ g DNA required

12 sample concurrently



# Original Paper

Future Medicine Ltd

Epigenomics

Volume 1, Issue 1, October 2009, Pages 177-200

<https://doi.org/10.2217/epi.09.14>

Technology Report

## Genome-wide DNA methylation profiling using Infinium<sup>®</sup> assay

Marina Bibikova<sup>1</sup>, Jennie Le<sup>1</sup>, Bret Barnes<sup>1</sup>, Shadi Saedinia-Melnyk<sup>1</sup>, Lixin Zhou<sup>2</sup>, Richard Shen<sup>1</sup> & Kevin L Gunderson<sup>1,†</sup>

<sup>1</sup>Illumina, Inc., 9885 Towne Centre Dr., San Diego, CA 92121, USA.

[kgunderson@illumina.com](mailto:kgunderson@illumina.com)


<sup>2</sup>State Key Laboratory of Medical Genetics, Central South University, 110 Xiangya Road, Changsha, Hunan 410078, PR China

<sup>††</sup> *Author for correspondence*

# Epigenomics




## Is it reliable ?




## Is it reliable ?

**A** *CD44* (four CpG sites in six samples)




**B** 113 CpG sites in Jurkat cell line



# Other comparisons

**Figure 1: Outline of the DNA methylation technology comparison.**

From: [Quantitative comparison of genome-wide DNA methylation mapping technologies](#)




Four methods for DNA methylation mapping were compared on two pairs of samples. The resulting 16 DNA methylation maps were bioinformatically analyzed and benchmarked against each other. In addition, clonal bisulfite sequencing was performed on selected genomic regions to validate DNA methylation differences that were detected exclusively by one method.

# Genome coverage

**Figure 4: Genomic coverage of MeDIP-seq, MethylCap-seq, RRBS and Infinium.**

From: [Quantitative comparison of genome-wide DNA methylation mapping technologies](#)




# Illumina 27K Array

## Promises

- Reliability
- Sample throughput

## Challenges




# Illumina 27K Array

## Promises

- Reliability
- Sample throughput

## Challenges

- Coverage (limited to 27k CpG sites; 0.1% of total CpGs)



\* 450K (1.7% of total CpGs), EPIC (3% of total CpGs)

# Illumina 27K Array

## Promises

- Reliability
- Sample throughput

## Challenges

- Coverage (limited to 27k CpG sites;  
0.1% of total CpGs)

Table 1: Markers of the HumanMethylation27BeadChip. <sup>3</sup>		
Type of target	CpG sites present	Avg # of CpG sites per target
RefSeq Genes	14,475	1.9 sites
Well-annotated genes described in the NCBI CCDS database (Genome build 36)	12,833	1.9 sites
Methylation hotspots in cancer genes	144	7.6 sites
Cancer-related targets	982	1.9 sites
miRNA promoters	110	2.3 sites



# Illumina 27K Array

## Promises

- Reliability
- Sample throughput

## Challenges


- Coverage (limited to 27k CpG sites; 0.1% of total CpGs)
- Sources of biases

Table 1


Potential sources of bias for different DNA profiling approaches.

Source of bias	Bisulfite conversion-based		Microarray-based	Endonuclease digestion-based	Affinity enrichment-based	
	WGBS	RRBS	Illumina DNA BeadChip	MRE-seq	MeDIP-seq	MBD-seq
Incomplete bisulfite conversion	✓	✓	✓			
Post-bisulfite conversion PCR	✓	✓	✓			
Cross-hybridization			✓			
Fragment size variation				✓		
Copy number variation					✓	✓
CpG density					✓	✓


## DNA methylation technology ...



## DNA methylation technology ...



## DNA methylation technology ...



## DNA methylation technology ...

⚠  
Modifications can be both  
stable and dynamic

⚠  
Data Integration is  
still crucial

DNA methylation studies

